

PREPARATION OF SSSSS DNA

Adapted from BV's Cookbook, TCH 1/27/02

- 1) Combine 3-3.3 g Salmon Sperm DNA (Sigma # S3126), 640 ml ddH₂O, and 40 ml 4M NaOH in a 1 liter screw cap glass flask.
- 2) Boil 20 minutes, shaking the flask occasionally. The DNA won't go into solution until it is boiled.
- 3) Put on ice until RT or colder.
- 4) Pour into a 2 liter flask. Add 75 ml 1M Tris pH 7.5 and 300 ul 0.5% Phenol Red.
- 5) pH to about 7.0 by adding about 14 ml concentrated HCl. Use phenol red as an indicator, since it turns orange/yellow at ~ pH 7. Also, check with pH paper.
- 6) Add 800 ml EtOH. Store overnight at 4°.
- 7) Pour into four 450 ml centrifuge tubes and spin 1 hour at 6K in GS3 rotor.
- 8) Pour off sup, then wash with 50 ml 70% EtOH per tube.
- 9) Dry inverted for about 30 minutes RT. Combine the pellets into one tube, then dry in a vacuum for about 1.5 hours.
- 10) Dissolve the pellet in 400 ml TE 10/2 pH 7.5 (20 ml 1 M Tris pH 7.5, 20 ml 0.2M EDTA, qs to 2.0 liters). Then transfer to a 2 liter flask and add 400 ml more TE.
- 11) Boil 20 minutes, then cool on ice.
- 12) Measure concentration by OD 260. First, dilute 1:100. Then, take OD 260 and OD 280. An OD 260 of 1.0 is about 43 ug/ml. OD 260/OD280 should be about 1.9.
- 13) Add TE (usually about 500 ml) to give a final concentration of 2 mg/ml.