

PROTOCOL FOR PC-8 EXTRACTION OF DNA SAMPLES

by TCH 1/26/02

1. Prepare your DNA solution as the follows:

| | |
|---------------------------|---|
| DNA solution: | 200ul (if <200ul, add ddH ₂ O to q.s. 200ul) |
| 7.5M NH ₄ Oac: | 100ul |
2. Add 250-300ul of phenol/chloroform, pH 8.0 (i.e., PC-8, available from Fisher or VWR). Vortex 15 sec and spin 1 min in microfuge (**Note:** PC-8 is in the bottom phase; always wear gloves).
3. Transfer the aqueous phase to a labeled fresh tube containing an equal volume of PC8. Vortex 15 seconds, microfuge 1-2 minutes (**Note:** It's best not to transfer all the aqueous phase, as this often gets the interphase junk mixed in; Usually 1-2 extractions are sufficient, but the cardinal rule is that there is no visible junk at the interphase after the final PC8 extraction).
4. Transfer the aqueous phase to a 1.7ml tube. Add 3ul seeDNA and 600-700ul 100% ethanol. Mix well and spin at top speed for 5 min in room temperature (**Note:** If you're purifying very short DNA, e.g., <100 bp, add three volumes, i.e., 900ul, of 100% ethanol).